





## DP 150

### ABOUT

The DP 150 has taken a sometimes arduous test protocol and reduced it down to its simplest elements. There is no longer any need to use color metric reference sheets or surfactants when developing a dye solution for the ASTM F-1929 protocol. Also, with the kit's included reference guide, you can quickly see what a seal failure looks like. In house packaging testing and inspection performed often and throughout the day is the new expectation under this heavily regulatory requirement for medical device packaging.

The DP 150 kit is affordable, so additional kits can be purchased throughout the facility to facilitate pouch and machine performance testing throughout the firm, including manufacturing, quality department, maintenance and even management.





## DP 150

### **DP 150 KITS CONTAINS THE FOLLOWING:**

- 1 - 75ml bottle of premixed testing dye
- 20 - disposable droppers
- Instructions for proper use and interpretation of results
- Material Data Safety sheet

### **DP 150 PRO KITS CONTAIN THE FOLLOWING:**

- 30 - premixed, individual foil dye pouches that can be placed inside of a pouch before sealing
- Instructions for proper use and interpretation of results
- Material Data Safety sheet