





## SPC 1 Specimen Cutter

### ABOUT

The SPC 1 peel test specimen cutting guide is ideal for any peel tensile test. Designed to work in concert with the PTT 100V peel testing unit, this specimen cutting guide boasts incredibly durable construction. The cutting guide will allow you to standardize the sample size used for test purposes that will deliver more consistent results, particularly when used with our specimen-loading tool. When performing seal strength test under the ASTM f-88 guidelines specimen size, and specimen repeatability is paramount for gathering exacting data.





## SPC 1 Specimen Cutter

### SPECIFICATIONS

- **Weight:** 3 pounds
- **Table Dimensions:** Width = 220mm, Length = 150mm, Height = 35mm
- **Sample Dimensions:** Width = 25mm, Length = 75mm
- Stainless steel base
- Heavy duty cutting mat
- Stainless steel cutting guide attached with hinges for easy sample placement and removal.
- Rubber feet prevent slipping during usage
- Alignment guide for seal position