



VIU

ABOUT & SPECIFICATIONS

This VIU allows the operator to perform a complete visual inspection of the seal area using our patented low angle of light and 3X magnification to better detect any flaws in the seal that might not present themselves to standardized seal strength testing.

- **Power Supply:** AC100-240V 50Hz/60Hz DC24V AC Adaptor
- **Power Consumption:** 1W
- **LED:** 3pc white LED
- **Lens:** Proprietary lens design allows 3X magnification
- **Data Output:** USB port
- **Photo Sensor:** Detects the presence of pouches for visual inspection
- **Operating Temperature Range:** 41°F – 104°F
- **Operating Humidity:** 30% - 80%
- **Dimensions:** 250mm wide x 195mm deep x 95mm high
- **Weight:** Approximately 5.5 pounds

